

SCIOLOGEX Syringe Tips

For Repeating Pipettes

The Perfect Tips for Your Repeating Pipette.

SCIOLOGEX syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- **Compatible with Most Pipettes:** SCIOLOGEX Syringe Tips are compatible with repeating pipettes using standard tips.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** SCIOLOGEX syringe tips are made entirely from high-quality virgin polyethylene and polypropylene. (Piston of 0.1 mL tip is liquid crystal polymer-LCP).
- **Available sterile:** SCIOLOGEX Sterile Syringe Tips are individually wrapped, and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.
- **Batch Certificate of Performance** which documents tip accuracy and precision testing is supplied with each pack of SCIOLOGEX Syringe Tips



Sterile, individually wrapped

Capacity,	Pack of	Cat. No.
0.5 mL	100	70230-
1.25 mL	100	70230/
2.5 mL	100	702389
5 mL	100	70*+91
12.5 mL	100	70*+93
25 mL (one adapter included)	25	70*+95
50 mL (one adapter included)	25	70*+97
25 mL and 50 mL Adapter	1	702398

Non-sterile, bulk packed

Capacity,	Pack of	Cat. No.
0.5 mL	100	702371
1.25 mL	100	702373
2.5 mL	100	702375
5 mL	100	702377
12.5 mL	100	702379
25 mL (one adapter included)	50	702381
50 mL (one adapter included)	25	702383
25 mL and 50 mL Adapter	1	702398

*The Eppendorf® Repeater® Plus, Repeater® M4, Repeater® Pro, Stream, and XStream are designed to preclude their operation with any tips other than Eppendorf® Combitips® advanced.

Toll Free: 877-724-5643 www.SCIOLOGEX.com

