



### **SCI-Filler**

# **Motorized Pipette** Filler

### **Features**

- Compatible with most of the plastic and glass pipettes from 0.1 -100mL
- Eight speeds of selection for aspiration and dispensing different liquids
- Large LED display showing low battery warning and speed settings
- Enables single-handed operation with minimum effort
- Light and ergonomic design providing easy usability
- High capacity Li-ion battery enable long runtime of operation
- Powerful pump fills a 25mL pipette in<5 seconds
- 0.45µm replaceable hydrophobic filter
- Rechargeable during use

Aspiration and dispensing operation buttons

Large LED display showing low battery warning and speed settings

Select the desired coarse speed by pressing the button +/-





Operating button for motor and gravity dispensing



Hydrophobic filter 0.45µm

Silicone adapter for 0.1-100mL

### **Specifications**

Aspiration Speed	25mL<7s (8 shift)
Dispensing Speed	Motor (8 shift)/Gravity
Battery	Lithium-ion
<b>Battery Service Life</b>	More than 8 Hours of Intermittence Use
Charging Time	2-3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45μm Hydrophobic
Weight	200g

## **Pipette Filler**

## **SCI-Filler2**

**Motorized Pipette Filler** 

SCI-Filler pipette filler has an excellent ergonomic comfortable and design, suitable for almost all commonly used glass or plastic serological pipettes.

#### **Features**

- Enables single-handed operation with minimum effort
- Compatible with most of the plastic and glass pipettesfrom 0.1-100mL
- Easily adjustable speed control
- 0.45µm replaceable hydrophobic filter
- Stable and reliable performance, elegant appearance
- High capacity Li-ion battery enable long operation time.







### **Specifications**

Aspiration Speed	25ml<7s
Dispensing Speed	Motor/Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2-3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45µm Hydrophobic
Weight	200g

## **Pipette Controllers & Pumps**

### **SCI-Pet**

## **Pipette Controller**

#### **Features**

- Ease of use
- Precise pipetting control
- Robust and light weight
- Compatible with most of the plastic and glass pipettes from 0.1 -100mL
- 3µm replaceable hydrophobic filter
- Easy cleaning and maintenance



Compress the suction bellows before attaching the pipette

Button to expel residual liquid when using blow-out pipette

Aspiration and dispensing comfortably achieved with a sensitive lever

### Levo E

## **Pipette Pump**

### **Features**

- Volume capacity 2mL, 10mL and 25mL
- Colour-coded by volume with green, blue and red
- Thumb wheel guarantees precise operation
- Resistant to acids and alkalis
- Easy for cleaning and maintenance

3µm replaceable Hydrophobic filter

Silicone adapter for 0.1-100mL
pipette tip, fully autoclavable

Designed for single-handed operation, Pipette Pump identifies itself with a thumbwheel, allowing precise aspiration or dispensing. The feature of side lever enables complete liquid dispensing, without a droplet remaining.

