



Designed for easy operation, Adam Equipment's MDW digital scales meet the needs of medical professionals and health practitioners worldwide. MDW digital scales feature an adjustable height rod to provide simultaneous height and weight readings, and some models also provide body mass index (BMI) measurements. All MDW digital scales feature a built-in battery and easy-to-use, color-coded keypad. To ensure stability and patient security, MDW digital scales feature rugged construction with a sturdy metal base and a non-slip weighing platform.

Affordable, easy-to-operate scale measures height and weight



Features

- BMI (body mass index) function available on the MDW 300L
- The MDW 300L offers RS-232 interface for fast communication with printers and computers
- Hold function retains the weight on the display
- Includes AC adapter and rechargeable battery
- Height measurements from 60cm to 212cm
- Integrated wheels permit easy movement (not available on MDW 250L)
- Anti-slip platform ensures safety when stepping on and off the scale
- Height rod facilitates measuring children and adults while weighing
- Metal construction for base and pillar

Applications

- Weighing

About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

MDW Digital Physician Scales

Model	MDW 250L	MDW 300L
Capacity	550lb / 250kg	660lb / 300kg
Readability	0.2lb / 0.1kg	0.1lb / 0.05kg
Repeatability (S.D.)	0.2lb / 0.1kg	0.2lb / 0.1kg
Linearity (+/-)	0.2lb / 0.1kg	0.2lb / 0.1kg
Pan Size	14.8"x10.8" / 375x275mm	
Weighing Units	kg, lb	
Stabilization Time (sec)	3	
Interface	N/A	RS-232
Body Mass Index [BMI]	No	Yes
Height Rod	Yes	
Calibration	External calibration	
Display	LED with 0.9" / 22mm -high digits	
Power Supply	Internal rechargeable battery / 12VDC 500mA adapter	
Operating Temperature	5° to 40°C / 41° to 104°F	
Overall Dim.	21.1"x10.8"x37" / 535x275x940mm (wxdxh)	22.4"x10.8"x51.6" / 570x275x1310mm (wxdxh)
Net Weight	33.1lb / 15kg	28.7lb / 13kg

Accessories

Item number	Description
3014011014	RS-232 cable M-F
3074010507	RS-232 to USB adapter (must order corresponding RS-232 cable)
1120011156	ATP thermal printer
3126011263	ATP thermal printer paper
3126011281	ATP thermal printer paper (pack of 10)
700660290	Calibration certificate

Key Accessories



ATP thermal printer
 Easy-to-use features combined with quick set-up equals a truly versatile thermal printer. A single push-button feed key and three easy-to-read LED indicators make operation simple.

Supplied by