

Sheldon Forced Air Oven Model SMO28-2

27.5 Cubic Ft. Upright Temperature Range: To 260°C

Specifications:

| | |
|----------------------|---|
| * Capacity | 27.5 cu. Ft. |
| * Temperature Range | 15°C above ambient to 260°C |
| * Dimensions | |
| * Exterior w x d x h | 37.5 x 35.0 x 78.3" 95.3 x 88.9 x 198.8 (cm) |
| * Interior w x d x h | 30.8 x 25.0 x 61.0" 78.2 x 63.5 x 157.5 (cm) |

Design Information:

The SMO28-2 is the largest general purpose oven on the market (also consider the SMO38HP-2 high performance ovens for the best uniformity and functionality). The independent overtemperature safety feature is user adjustable and provides added safety. The turbo blower-assisted forced airflow design coupled with a heavy-duty blower motor offers excellent temperature uniformity and fast recovery times. Airflow is directed evenly across all shelves and samples to ensure consistent drying, curing and baking results.

Features:

| | |
|----------------------------|-------------------------|
| Stainless Steel Shelves | Yes |
| 3" Adjustable Exhaust Port | Yes |
| No Exposed Heating | Yes |
| Wrap Around Insulation | Yes |
| Digital Timer | Yes |
| Temperature Uniformity | 3.0°C |
| Temperature Stability | 0.2°C |
| Shelving | 6 Supplied (16 Maximum) |
| Voltage | 230/60/1 |
| Amperage | 20 |
| Approval | CE Approved |
| Shipping Weight | 495 lbs. |
| Warranty | Two Year |

Specifications subject to change without notice

