

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-10SD

10 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Solid Door

Specifications:

* Capacity	10 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	24" x 22" x 60" 60 x 55 x 152 (cm)
* Interior w x d x h	19.5" x 17" x 50" 49 x 43 x 127 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Door	Single solid, locking
Access Port	1/2" Diameter
Interior Finish	White coated steel
Exterior Finish	White coated steel
Legs	Leveling Legs
Shelves	Five Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Compressor	1/5 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	1.8 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	300 lbs.
Warranty-(Domestic)	Two year parts and labor.

