

## MODEL DHK4-14SD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>14 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>395 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>35°F to 46°F</b>

Dimensions Width x Depth x Height

Exterior [ Inches ] 25.75" x 26.5" x 78.5"

Interior [ Inches ] 23" x 21" x 53.25"

Exterior [ Centimeters ] 65 x 67 x 199

Interior [ Centimeters ] 58 x 53 x 135

Compressor Hermetic Compressor

Air Circulation Forced Air

Defrost Method Automatic

Refrigerant R600a

Exterior Steel

Interior Polystyrene ( HIPS )

Insulation US EPA and SNAP approved

Lockable Door Yes, Keyed

Porthole Yes, 1" Diameter

Casters 4, with Leveling Feet

Storage 6 Adjustable Shelves, 1 Basket

Shipping Weight 250 lbs.

#### Electrical

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage Line | Running 15 Amp Dedicated | 3.77A

Power Consumption 3.71 kW/d

NEMA Configuration NEMA 5-15

Certification(s) UL Listed

#### Controller

Display Type Microprocessor, Digital

Alarms High/Low Temp, Door Ajar

Power Failure, Low Battery

Sensor Failure, USB Failure

Yes, Adjustable Intervals

Data Logging Yes, via USB, PDF Format


Data Download Yes, Display and Reset

Min / Max Temperature Yes, Dry Contacts

Alarm Relay Yes, Only Powers Alarm

Backup Battery

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor