

SO-LOW EXPLOSION PROOF FREEZER MODEL DHH20-17SDFX MANUAL DEFROST

Meets NFPA and OSHA guidelines for Hazardous Locations, Class I, Division II, Group C & D.

17 Cubic Ft. Upright -25°C / Solid Door

Specifications:

* Capacity	17 cubic ft. / 481 liters
* Temperature	-15°C to -25°C
* Dimensions	
* Exterior w x d x h	30.75" x 28.50" x 67.75"
* Interior w x d x h	26.25" x 22.50" x 58.25"

Features:

*Thermostat	*Microprocessor temperature controller *Remote alarm contacts *Audible and visual high/low temperature alarms *Alarm validation (audible, visual, and remote alarm)
*Compressor	1/2 Hp Compressor
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	6 fixed shelves
*Insulation:	High density urethane foam
*Refrigerant:	R-290A
*Air Circulation	Cold wall
*Defrost	Manual
*Voltage	115-60-1
*Running Amps	3.5
*Supply Plug	Not supplied (hardwire)
*Shipping weight	220 lbs.
*Approvals	ETL Listed
*Warranty-(Domestic)	Two year parts and labor.



Specifications subject to change without notice