

SO-LOW EXPLOSION PROOF FREEZER MODEL DHH20-20SDFX MANUAL DEFROST

Meets NFPA and OSHA guidelines for Hazardous Locations, Class I, Division II, Group C & D.

20 Cubic Ft. Upright -25°C / Solid Door

Specifications:

* Capacity	20 cubic ft. / 566 liters
* Temperature	-15°C to -25°C
* Dimensions	
* Exterior w x d x h	30.75" x 29" x 74.50"
* Interior w x d x h	25.75" x 21.75" x 65"

Features:

*Thermostat	*Microprocessor temperature controller
	*Remote alarm contacts
	*Audible and visual high/low temperature alarms
	*Alarm validation (audible, visual, and remote alarm)
*Compressor	1/2 Hp Compressor
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	6 fixed shelves
*Insulation:	High density urethane foam
*Refrigerant:	R-290A
*Air Circulation	Cold wall
*Defrost	Manual
*Voltage	115-60-1
*Running Amps	6
*Supply Plug	Not supplied (hardwire)
*Shipping weight	240 lbs.
*Approvals	ETL Listed
*Warranty-(Domestic)	Two year parts and labor.



Specifications subject to change without notice