

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR/FREEZER*

### *MODEL DHH-9RFFMS*

*Units meet NFPA guidelines 45 and 70*

**Refrigerator 4.6 cu. ft. (130 liters)-Freezer 4.4 cu. ft. (125 liters)**

**Upright – Temperature Range: 1°C to 10°C (R) / -15°C to -25°C (F) Two Doors**

#### **Specifications:**

* Capacity	9 cubic ft. / 255 liters
* Temperature Range	1°C to 10°C / -15°C to -25°C
* Dimensions	
* Exterior w x d x h	23 5/8" x 24 1/4" x 65 3/8" 60 x 61.60 x 168.61 (cm)
*Interior w x d x h	20.25" x 17" x 27" - Refrigerator 51.44 x 43.18 x 68.58 (cm) 19" x 16" x 26" - Freezer 48.26 x 40.64 x 66.04 (cm)

#### **Features:**

Temperature Control	Two - Digital display
Temperature Alarm	Hi/Low temperature alarm With audible, and visual alarm Alarm relay dry contacts
Doors	Two solid, locking
Legs	Four leveling feet
Shelves	Three adjustable / Refrigerator Two Fixed/ Freezer
Stacking	Stacking Kit
Insulation	Urethane foam
Refrigerant:	R600A CFC free
Compressor	Two
Thermostat	Adjustable
Refrigerator Defrost	Cycle, automatic
Freezer Defrost	Manual
Voltage: Domestic	115/60/1
Supply Plug	(2)5-15P NEMA
Shipping Weight	205 lbs.
Warranty-(Domestic)	Two year parts and labor Five year compressor warranty



Specifications subject to change without notice