

SO-LOW FLAMMABLE MATERIAL STORAGE MANUAL DEFROST REFRIGERATOR/FREEZER MODEL DHH-30RFFMS

Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas

Refrigerator 25 cu. ft. (707 liters)-Freezer 5.0 cu. ft. (141 liters)

Upright – Temperature Range: 1°C to 10°C (R) / -15°C to -25°C (F)

Specifications:

* Capacity	30 cubic ft. / 849 liters
* Temperature	1°C to 10°C / -15°C to -25°C
* Dimensions	
* Exterior w x d x h	35 ¾" x 35 ½" x 72 ½" 90.81 x 90.17 x 184.15 (cm)
* Interior w x d x h	29" x 22.25" x 12" - Freezer 73.66 x 56.52 x 30.48 (cm) 32" x 22.5" x 50.25" - Refrigerator 81.28 x 57.15 x 127.64 (cm)

Features:

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	3 adjustable – Fridge only
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-404A
*Compressor	1/3 hp
*Air Circulation	Gravity
*Evaporator	Hot wall
*Defrost	Manual
*Voltage	115/60/1
*Running Amps	6
*Supply Plug	5-15 NEMA
*Shipping weight	396 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.

