

## *SO-LOW LABORATORY FREEZER TO -20°C*

### *MODEL DHW20-30*

**30 Cubic Ft. (849liters) Upright**

**Temperature Range: -20°C**

**Solid Door**

#### **Specifications:**

* Capacity		30 cubic ft. / 849 liters
* Temperature Range		-20°C
* Dimensions		
* Exterior	w x d x h	35.5" x 35" x 72"
* Internal	w x d x h	32" x 22.5" x 65"

#### **Features:**

- Temperature Control
- Manual Defrost
- R404A Refrigeration System
- 4 Fixed Shelves
- Right Hand Swing Door
- Heavy Duty Lift off Hinges
- Stainless Steel Handle
- Key Lock
- 115 Volt, 60 Hertz, 1 Phase
- One Year Warranty

