

SO-LOW LABORATORY FREEZER TO -20°C

MODEL DHW20-30MDP

30 Cubic Ft. (849liters) Upright

Temperature Range: -20°C

Solid Door

Specifications:

* Capacity		30 cubic ft. / 849 liters
* Temperature Range		-10°C to -20°C
* Dimensions		
* Exterior	w x d x h	35.75" x 34.25" x 72.5"
* Internal	w x d x h	32" x 22.5" x 65" 81.28 x 57.15 x 165.1

Features:

CONTROL

Digital temperature control

ALARM

Alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting

REMOTE ALARM RELAY

A relay for a remote alarm hook-up is also provided.

Manual Defrost

R404A Refrigeration System

Heavy Duty Lift Off Hinges

4 Fixed Shelves

Right Hand Swing Door

Stainless Steel Handle

Key Lock

115 Volt, 60 Hertz, 1 Phase

6 Amps

One Year Warranty

