

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB AND PHARMACY REFRIGERATOR FREEZER COMBINATION MODEL DH-20RF

Refrigerator 9.1 cu. ft. -Freezer 9.1 cu. ft.

Upright – Temperature Range: 8°C (R) / -20°C (F) Two Doors

Specifications:

* Total Capacity	18.2 cubic ft.
* Temperature (Fridge)	2°C to 8°C
* Temperature (Freezer)	0°C to -20°C
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	27.5" x 34" x 81.5" 70 x 86 x 207 (cm)
* Fridge Interior w x d x h	28 x 28 x 20"
* Freezer Interior w x d x h	28 x 28 x 20"

Features:

Control	Digital Display with Hi/Low Alarm Alarm Relay Dry Contacts For each compartment
Defrost	Auto
Doors	Two Solid, Locking
Casters	Four, Lockable
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Shelves	Epoxy Coated / 2 Per Compartment
Insulation:	Polyurethane
Evaporator	Fin and tube
Compressors	1/4 hp (fridge), 1/3 hp (freezer)
Voltage:	115/60/1
Running Amps	4.5 (fridge), 7.5(freezer)
Supply Plugs	(2) NEMA 5-15P UL listed
Shipping Weight	550 lbs.
Warranty	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice