

## ***SO-LOW LAB AND PHARMACY REFRIGERATOR FREEZER COMBINATION MODEL DH-48RF***

**Refrigerator 24 cu. ft. -Freezer 24 cu. ft.**

**Upright – Temperature Range: 8°C (R) / -25°C (F) Two Doors**

### **Specifications:**

\* Total Capacity 48 cubic ft.  
\* Temperature (Fridge) 2°C to 8°C  
\* Temperature (Freezer) 0°C to -25°C

#### **\* Dimensions**

\* Exterior w x d x h 52 x 34.75 x 81.5”  
132 x 88 x 207”

\* Fridge Interior w x d x h 22 x 28 x 60”

\* Freezer Interior w x d x h 22 x 28 x 60”

### **Features:**

Control	Digital Display with Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory For each compartment
Defrost	Auto
Doors	Two Solid, Locking
Casters	Four, Lockable
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Shelves	Epoxy Coated / Three Per Compartment
Insulation:	Polyurethane
Evaporator	Fin and tube
Compressor	1/4 hp - 1/2 hp
Voltage:	115/60/1
Running Amps/Breaker	5.7/15 – 6.0/15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	625 lbs.
Warranty-(Domestic)	3 year parts and labor Compressor additional 4 years, no labor.

