

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-23GD

23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass Door

Specifications:

| | |
|----------------------|---|
| * Capacity | 23 cubic ft. |
| * Temperature Range | 2°C to 8°C / 36°F to 46°F |
| * Dimensions | |
| * Exterior w x d x h | 27" x 32" x 81" 68 x 81 x 206 (cm) |
| * Interior w x d x h | 24" x 28.5" x 55" 61 x 72 x 140 (cm) |

Features:

| | |
|----------------------|--|
| Temperature Control | Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory |
| Access Port | 1/2" Diameter Top Back Wall |
| Doors | Single glass, locking |
| Interior finish | White powder coat |
| Casters | Four, locking |
| Shelves | Five |
| Interior Lighting | LED light |
| Insulation | Polyurethane foam |
| Refrigerant: | R290 CFC free |
| Compressor | 1/3 hp |
| Air Circulation | Forced air |
| Defrost | Cycle Defrost |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 3.0/15 |
| Supply Plug | 5-15P NEMA UL listed |
| Shipping Weight | 440 lbs. |
| Warranty-(Domestic) | 3 year parts and labor Compressors additional 4 years, no labor. |

