

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-38SGD

38 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Two Glass Sliding Doors

Specifications:

* Capacity	38 cubic ft. / 1076 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	43.5" x 32" x 81" 110 x 81 x 205 (cm)
* Interior w x d x h	40.5 x 27" x 55" 102 x 68 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Double glass – sliding Lockable
Casters	Four, locking
Shelves	Ten, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps	5.0/15
Supply Plug	5-15 NEMA UL listed
Shipping Weight	550 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

