

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-66SGD

66 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Three Sliding Glass Doors

Specifications:

* Capacity	66 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	75" x 32" x 81" 190 x 81 x 206 (cm)
* Interior w x d x h	72" x 27" x 55" 182 x 69 x 140 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Three sliding glass, locking
Casters	Four, Locking
Shelves	Fifteen, Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/2 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.2 / 15
Supply Plug	5-15P NEMA
Shipping Weight	850 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice