

SO-LOW LABORATORY & PHARMACY PASS-THRU REFRIGERATOR MODEL DH4-25PT

25 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass doors on front and back

Specifications:

* Capacity	25 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	26" x 37.5" x 81.5" 66 x 95 x 207 (cm)
* Interior w x d x h	23" x 30" x 59" 58 x 76 x 149 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm Alarm Relay Dry Contacts Min/Max Memory
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Door	Glass doors front and back, lockable Solid door option available
Casters	Four, locking
Shelves	Three, Adjustable
Interior Light	One
Insulation:	Polyurethane foam
Refrigerant:	R-290
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage:	115/60/1
Running Amps	4.2
Supply Plug	5-15P NEMA UL listed
Shipping Weight	500 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice.