

SO-LOW LAB & PHARMACY REFRIGERATOR

MODEL DH4-80PT

80 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass doors on front and back

Specifications:

* Capacity	80 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	78" x 37.5" x 81.5" 198 x 95 x 207 (cm)
* Interior w x d x h	75" x 30" x 59" 190 x 76 x 149 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm Alarm Relay Dry Contacts Min/Max Memory
Construction	Stainless Interior and Exterior
Door	Glass doors front and back, lockable Solid door option available
Casters	Four 3" swivel, locking
Shelves	Nine, Adjustable
Interior Light	Two Fluorescents
Insulation:	Polyurethane foam
Refrigerant:	R-290
Compressor	1 hp
Air Circulation	Forced air
Evaporator	Fin and tube
Defrost	Cycle Defrost
Voltage:	115/60/1
Running Amps	14.0
Supply Plug	5-15P NEMA
Shipping Weight	850 lbs.
Warranty-(Domestic)	3 year parts and labor. Compressors additional 4 years, no labor.



Specifications subject to change without notice.