

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-23GDSS*

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Glass Door**

### **Specifications:**

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	27.25" x 33.75" x 81" 69 x 85 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

### **Features:**

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory
Doors	Single glass, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four, locking
Shelves	Three
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice