

MODEL MV4-2UCRGDDA

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

| | |
|---------------------------------|---------------------|
| Capacity [Cu. Ft.] | 2.64 Cu. Ft. |
| Capacity [Liters] | 75 Liters |
| Temperature Range [°C] | 2°C to 8°C |
| Temperature Range [°F] | 35°F to 46°F |

| | |
|--------------------------|------------------------|
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 21.25" x 22" x 30" |
| Interior [Inches] | 17.5" x 17.5" x 21.75" |
| Exterior [Centimeters] | 54 x 56 x 76 |
| Interior [Centimeters] | 44 x 44 x 54 |

| | |
|-----------------|--------------------------------|
| Compressor | Hermetic Compressor |
| Air Circulation | Forced Air |
| Defrost Method | Automatic |
| Refrigerant | R600a |
| Exterior | Steel |
| Interior | Polystyrene (HIPS) |
| Insulation | US EPA and SNAP approved |
| Lockable Door | Yes, Keyed |
| Porthole | Yes, 1" Diameter |
| Casters | 4, with Leveling Feet |
| Storage | 4 Adjustable Shelves, 1 Basket |
| Shipping Weight | 150 lbs. |


Electrical

| | |
|-------------------------|----------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 1.24A |
| Power Consumption | 1.24 kW/d |
| NEMA Configuration | NEMA 5-15 |
| Certification(s) | UL Listed |

Controller

| | |
|-----------------------|---|
| Display Type | Microprocessor, Digital |
| Alarms | High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure |
| Data Logging | Yes, Adjustable Intervals |
| Data Download | Yes, via USB, PDF Format |
| Min / Max Temperature | Yes, Display and Reset |
| Alarm Relay | Yes, Dry Contacts |
| Backup Battery | Yes, Only Powers Alarm |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor