

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW UNDERCOUNTER FREEZER

MODEL MV20-1.5UCF

1.5 Cubic Ft. (42 liters) Upright

Temperature Range to -20°C

Solid Door

Specifications:

* Capacity	1.5 cubic ft. / 42 liters
* Temperature Range	-20°C
* Dimensions	
* Exterior w x d x h	14-7/8" x 22-1/8" x 24 1/4" 38 x 56 x 62 (cm)
* Interior w x d x h	11-1/4" x 17-5/8" x 12" 29 x 45 x 30 (cm)

Features:

COMMERCIAL U.L. LISTED

MANUAL DEFROST CYCLE

FORCED AIR VENTILATION
THROUGH FRONT GRILLE

ADJUSTABLE THERMOSTAT CONTROL

NO ADDITIONAL VENTILATION
CLEARANCE REQUIRED TO BUILD
IN UNDERCOUNTER

404A REFRIGERATION SYSTEM

FLAT DOOR LINER

115 VOLT, 60 HERTZ, 1 PHASE

3.0 AMPS



Specifications subject to change without notice