

SO-LOW EXPLOSION PROOF REFRIGERATOR MODEL DHH4-17SDRX

Meets NFPA and OSHA guidelines for Hazardous Locations, Class 1, Division II, Group C & D.

17 Cubic Ft. Upright +4°C / Solid Door

Specifications:

| | |
|----------------------|---------------------------|
| * Capacity | 17 cubic ft. / 481 liters |
| * Temperature | 1°C to 10°C |
| * Dimensions | |
| * Exterior w x d x h | 30.75" x 28.5" x 67.75" |
| * Interior w x d x h | 26.25" x 20.75" x 58.50" |

Features:

| | |
|----------------------|---|
| *Thermostat | *Microprocessor temperature controller |
| | *Remote alarm contacts |
| | *Audible and visual high/low temperature alarms |
| | *Alarm validation (audible, visual, and remote alarm) |
| *Compressor | 1/4 Hp Compressor |
| *Doors | Single solid – key lock |
| *Legs | Four leveling |
| *Shelves | 5 shelves |
| *Insulation: | High density urethane foam |
| *Refrigerant: | R-600A |
| *Air Circulation | Fan |
| *Defrost | Cycle |
| *Voltage | 115-60-1 |
| *Running Amps | 3 |
| *Supply Plug | Not supplied (hardwire) |
| *Shipping weight | 240 lbs. |
| *Approvals | ETL Listed |
| *Warranty-(Domestic) | Two year parts and labor. |



Specifications subject to change without notice