

Containers

Product code	Basket model	Description	Dimensions	Suitable models
H060303020	-	Mesh basket	ø254 x H240mm	SK101C/111C
H060303038	-	Mesh basket	ø254 x H331mm	SK200C/210C
H060303037	-	Mesh basket	ø254 x H409mm	SK300C/310C
H060601028	-	Mesh basket	ø424 x H200mm	SQL810C/1010C
A990201303	-	Stainless rack	ø277 x H420mm	SK301C/311C
A990201304	-	Stainless rack	ø277 x H330mm	SK201C/211C
A990201305	-	Stainless rack	ø277 x H260mm	SK102C/112C
241085	OSQ-10	Mesh basket (pitch 2.5 mm)	ø190 x H159mm	SM201/210
241086	OSQ-20	Mesh basket (pitch 2.5 mm)	ø250 x H201mm	SM301/311/501/511
241087	OSM-60	Mesh basket (pitch 8.5 mm)	ø205 x H204mm	SM201/210, SN200C/210C
241088	OSM-70	Mesh basket (pitch 8.5 mm)	ø262 x H204mm	SM301/311/501/511, SN300C/310C/500C/510C
241089	OSM-80	Mesh basket (pitch 8.5 mm)	ø262 x H315mm	SM501/511, SN500C/510C
241090	OSR-40	Mesh basket with 2 stacking fittings	ø326 x H165mm	SQ500C/510C/810C, SM520/530/820/830
241091	OSQ-50	Mesh basket with 3 stacking fittings	ø246 x H162mm	SM501/511, SN500C/510C
241092	OSQ-40	Mesh basket with 2 stacking fittings	ø246 x H162mm	SM301/311, SN300C/310C
241093	OSQ-30	Mesh basket with 2 stacking fittings	ø168 x H162mm	SN200C/210C, SM201/210
241094	OSQ-80	Mesh basket with 2 perforated plates	ø200 x H590mm	SM501/511, SN500C/510C
241095	OSQ-70	Mesh basket with 1 perforated plate	ø260 x H390mm	SM301/311, SN300C/310C
241096	OSQ-60	Mesh basket with 1 perforated plate	ø200 x H390mm	SM201/210, SN200C/210C
241097	OSR-50	Mesh basket with 1 perforated plate	ø340 x H430mm	SQ500C/510C/810C, SM520/530/820/830
241099	OSQ-90	Mesh basket	ø332 x H196mm	SQ500C/510C/810C, SM520/530/820/830
241083	OSR-10	Stainless solid basket	ø210 x H200mm	SM201/210, SN200C/210C
241084	OSR-20	Stainless solid basket	ø270 x H200mm	SM301/311/501/511, SN300C/310C/500C/510C
241150	OSN-10	Stainless bucket	ø210 x H210mm	SM201/210, SN200C/210C
241151	OSN-12	Stainless bucket	ø270 x H210mm	SM301/311/501/511, SN300C/310C/500C/510C
241152	OSN-14	Stainless bucket	ø340 x H210mm	SQ500C/510C/810C, SM520/530/820/830

Mesh basket



Stainless bucket



Stainless solid basket




Plates and Inserts

Product code	Description	Suitable models
Q110603004	Stackable plate	SK200C/210C
Q110603005	Stackable plate	SK300C/310C
Q110603008	Stackable insert	SK200C/210C
Q110603009	Stackable insert	SK300C/310C




Stackable plate




Stackable insert

Shelves


Product code	Punching shape	Suitable sterilizer models
212095	Round punch shelf & bracket set	SI401/402
212246	Round punch shelf & bracket set	SK401
212266	Round punch shelf & bracket set	SI601/602, SK601/801/811



212095



212246




212266

Stands


Product code	Stand models	Suitable sterilizer models
212802	ONS60	SI401/402/601/602
211856	ON61	SK401/601
212348	OT42	SK401
212349	OT62	SK601



ONS60



ON61



OT42/62

