

SINCE 1889



LARGER CAPACITY LABORATORY USE STERILIZER

SQ500C/510C/810C

The large capacity and ergonomically designed SQ sterilizer makes it user friendly for sample handling

Features

- Comfortable accessibility and easy operation
- Easy to clean and maintain
- Mobile on wheels
- With high performance and safety functions
- Programmable or preset programs for commonly used sample types
- Reliable high pressure steam sterilization up to 135°C
- Cooling fan to shorten cool down time
- Various add-on options and wide selection of baskets/containers available



Easily manage the samples



Control panel

Specifications

Model	SQ500C/510C	SQ810C
Operating temp. range	45~135°C	
Max. operational pressure	0.255 MPa	
Cooling	Equipped with cooling fan	
Safety devices	Self diagnostic function, locking lid mechanism, safety pressure valve	
Internal dimensions	Dia 370 x 470 mm	Dia 370 x 750 mm
Internal capacity	50L	80L
Power source AC100-120V 50/60Hz Power source AC200-240V 50/60Hz	SQ500C 24.5~20.5A SQ510C 12.5~10.5A	SQ810C 12.5~10.5A
Weight	~95kg	~105kg