

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW PLATINUM SERIES REFRIGERATOR*

### *MODEL DHP4-23GD*

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Glass Door**

#### **Specifications:**

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	27" x 32" x 81" 68 x 81 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

#### **Features:**

Temperature Control	Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access
Access Port	1/2" Diameter Top Back Wall
Doors	Single glass, locking
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four, locking
Shelves	Five, Adjustable
Interior Lighting	One
Insulation	Polyurethane foam
Refrigerant:	R290A CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0 /15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty (Domestic)	3 year parts and labor. Compressors additional 4 years, no labor.



Specifications subject to change without notice