

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHP4-23SD*

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Solid Door**

### **Specifications:**

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	27" x 34" x 81" 68 x 86 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

### **Features:**

Temperature Control	Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access
Access Port	1/2" Diameter Top Back Wall
Doors	Single Solid, locking
Interior finish	White powder coat
Casters	Four, lockable
Shelves	Five
Insulation	Polyurethane foam
Refrigerant:	R290A CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty (Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

