

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-38SGD

38 Cubic Ft. Upright – Temperature Range: 2°C to 8°C

Two Glass Sliding Doors

Specifications:

* Capacity	38 cubic ft. / 1076 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	43.5" x 32" x 81" 110 x 81 x 205 (cm)
* Interior w x d x h	40.5 x 27" x 55" 102 x 68 x 140 (cm)

Features:

Temperature Control	Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access
Access Port	1/2" Diameter Top Back Wall
Interior	White powder-coated Steel
Doors	Double glass – sliding Positive seal, self closing, Lockable
Casters	Four, locking
Shelves	Ten
Interior Lighting	One fluorescent light
Insulation	Polyurethane
Refrigerant:	R290A
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps	5.0 /15 Amp circuit required
Supply Plug	5-15 NEMA
Shipping Weight	550 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

