

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW PLATINUM SERIES REFRIGERATOR MODEL DHP4-74SD*

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C  
Three Hinged Solid Doors**

### **Specifications:**

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	75" x 34" x 81" 190 x 86 x 206 (cm)
* Interior w x d x h	72" x 27" x 55" 183 x 69 x 140 (cm)

### **Features:**

Temperature Control	Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access
Access Port	1/2" Diameter Top Back Wall
Interior Doors	White powder-coated Three solid – hinged – key lock Positive seal, self closing
Shelves	Fifteen, Adjustable
Casters	Four, Locking
Interior Lighting	Two fluorescent lights
Insulation	Polyurethane CFC Free
Refrigerant	R290A CFC free
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	7.2/15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	850 lbs.
Warranty-(Domestic)	3 year parts and labor Compressor additional 4 years, no labor.



Specifications subject to change without notice