

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

***SO-LOW LAB & PHARMACY REFRIGERATOR  
 MODEL DHP4-45SGD***

**45 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Two Siding Glass Doors**

**Specifications:**

* Capacity	45 cubic ft. / 1286 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	52" x 32" x 81" 132 x 81 x 206 (cm)
* Interior w x d x h	49" x 27" x 55" 125 x 69 x 140 (cm)

**Features:**

Temperature Control	Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access
Access Port	1/2" Diameter Top Back Wall
Doors	Two sliding glass, locking
Casters	Four, Locking
Shelves	Ten
Interior Lighting	One fluorescent
Insulation	polyurethane foam
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0 /15
Supply Plug	5-15P NEMA
Shipping Weight	610 lbs.
Warranty-(Domestic)	3 year parts and labor Compressors additional 4 years, no labor.



Specifications subject to change without notice