

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHP4-49GD

49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C

Two Swinging Glass Doors

Specifications:

| | |
|----------------------|--|
| * Capacity | 49 cubic ft. / 1387 liters |
| * Temperature Range | 2°C to 8°C / 36°F to 46°F |
| * Dimensions | |
| * Exterior w x d x h | 52" x 32" x 81" 132 x 81 x 206 (cm) |
| * Interior w x d x h | 47.5" x 28.5" x 55" 120 x 72 x 140 (cm) |

Features:

| | |
|----------------------|---|
| Temperature Control | Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access |
| Access Port | 1/2" Diameter Top Back Wall |
| Doors | Two hinged glass, locking |
| Casters | Four, Locking |
| Shelves | Ten |
| Interior Lighting | One fluorescent |
| Insulation | Polyurethane foam |
| Refrigerant: | R290A |
| Compressor | 1/3 hp |
| Condenser | Fin and tube |
| Air Circulation | Forced air |
| Evaporator | Fin and tube, forced air |
| Defrost | Cycle Defrost |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 5.0 /15 |
| Supply Plug | 5-15P NEMA |
| Shipping Weight | 610 lbs. |
| Warranty (Domestic) | 3 year parts and labor Compressors additional 4 years, no labor. |

