

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHP4-49SD

49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C

Two Swinging Solid Doors

Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	52" x 34" x 81" 132 x 86 x 206 (cm)
* Interior w x d x h	47.5" x 28.5" x 55" 120 x 72 x 140 (cm)

Features:

Temperature Control	Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access
Access Port	1/2" Diameter Top Back Wall
Doors	Two hinged solid, locking
Casters	Four, Locking
Shelves	Ten
Interior Lighting	One fluorescent
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/3 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA
Shipping Weight	610 lbs.
Warranty (Domestic)	3 year parts and labor Compressors additional 4 years, no labor.

