

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHP4-66SGD

66 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Three Sliding Glass Doors

Specifications:

| | |
|----------------------|--|
| * Capacity | 66 cubic ft. |
| * Temperature Range | 2°C to 8°C / 36°F to 46°F |
| * Dimensions | |
| * Exterior w x d x h | 75" x 32" x 81" 190 x 81 x 206 (cm) |
| * Interior w x d x h | 72" x 27" x 55" 182 x 69 x 140 (cm) |

Features:

| | |
|----------------------|---|
| Temperature Control | Touch Screen Digital Display Hi/Low Alarm W/ Email and SMS notification Alarm Relay Dry Contacts Viewable Temperature graph Data logging downloadable via USB or FTP Multi-Level Security VNC Remote Access |
| Access Port | 1/2" Diameter Top Back Wall |
| Interior Finish | White coated steel |
| Exterior Finish | White coated steel |
| Doors | Three sliding glass, locking |
| Casters | Four, Locking |
| Shelves | Fifteen, Adjustable |
| Interior Lighting | LED Light |
| Insulation | Polyurethane foam |
| Refrigerant: | R290 CFC free |
| Compressor | 1/2 hp |
| Air Circulation | Forced air |
| Defrost | Cycle Defrost |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 6.2 / 15 |
| Supply Plug | 5-15P NEMA |
| Shipping Weight | 850 lbs. |
| Warranty-(Domestic) | 3 year parts and labor Compressors additional 4 years, no labor. |



Specifications subject to change without notice