

SO-LOW UPRIGHT FREEZER TO -40°C

MODEL NU40-25

25 Cubic Ft. (710 liters) Upright – Temperature Range: 0°C to -40°C

Upright Freezer

Specifications:

* Capacity	25 cubic ft. / 710 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior W x D x H	50" x 37" x 79.5" 127 x 93 x 201 (cm)
* Interior W x D x H	34" x 25" x 51.25" 86 x 63 x 130 (cm)

COLOR TOUCH SCREEN TEMPERATURE CONTROL

Touch Screen control

Displays set point and chamber temperature

Hi/Low Alarm W/ Email and SMS notification

Alarm Relay Dry Contacts

Power Failure Dry Contacts

Battery Back-up

Viewable Temperature graph

Data logging downloadable via USB or FTP

Multi-Level Security

VNC Remote Access

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible & visual signal when there is a mechanical or electrical failure. Relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low single stage system with one hermetic compressor. CFC & HCFC free refrigerants.

CONSTRUCTION

Chamber is 14-ga. Zinc coated galvanized steel.

Exterior is 16-ga. steel. Power coated cool gray.

VOLTAGE

20 AMP dedicated circuit required.

ETL electrical approval.

Specifications subject to change without notice

