



Model:

MS-1

Typology:

MEASURING MONOCULAR MICROSCOPE

Description:

Bench microscope for measuring and surface inspection.

Illumination	Penlight unit with xenon lamp. Powered by batteries (2x AAA 1,5V)
Focusing	By means of a rack and pinion controlled by a pair of knobs placed on both sides of the microscope.
Stand	Compound microscope 230 mm height
Eyepiece	WF10x/15.5 with field number 15.5. Micrometric with 0,1 mm division.
Objective	Achromatic 10X