


# MY-ISO5 Vertical Laminar Flow Cabinet - 110V

# MYSTAIRE®



Made in the USA



## GENERAL SPECIFICATIONS

<b>Construction Material</b>	<b>Work Surface</b>	Polypropylene
	<b>Shell</b>	Polycarbonate
	<b>Front Sash</b>	Polycarbonate
	<b>Filtration Housing</b>	ABS
	<b>Pipette Shelf</b>	Polycarbonate
<b>Air Quality</b>		Class 100 / ISO 5
<b>Controls</b>		ON/OFF Switch
<b>Blowers</b>		Brushless-Sparkless
<b>Lighting</b>		Fluorescent
<b>IV Bar</b>		Optional, Stainless Steel
<b>HEPA Filter Monitoring Package</b>		Audible / Visual Alarm

## DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-VLF24	24"	24"	30"	23"	23"	18"	21"	10"	20"
MY-VLF32	32"	24"	30"	31"	23"	18"	29"	10"	20"
MY-VLF32T	32"	24"	39"	31"	23"	27"	29"	11"	25"
MY-VLF48	48"	24"	30"	47"	23"	18"	45"	10"	20"
MY-VLF48T	48"	24"	42"	47"	23"	29"	45"	11"	25"

## TECHNICAL SPECIFICATIONS

Product #	Electrical		Volume of Filtered Air	Prefilters Required	Main Filters Required
	Voltage	Power			
MY-VLF24	110V	120 W	1	1	55 lbs.
MY-VLF32	110V	120 W	1	1	70 lbs.
MY-VLF32T	110V	120 W	1	1	80 lbs.
MY-VLF48	110V	220 W	2	2	95 lbs.
MY-VLF48T	110V	220 W	2	2	105 lbs.

Unit weight is an approximation\*

[www.mystaire.com](http://www.mystaire.com)

Phone: +1 (919) 229-8511 • Toll Free: (877) 328-3912 • Email: [contact@mystaire.com](mailto:contact@mystaire.com)