

Unix Spring Loaded Single Cartridge Holder (30 PSI Max)

UNIX-00

UNIX SPRING LOADED SINGLE CARTRIDGE HOLDER (30 PSI MAX)

Unix Spring Loaded Single Cartridge Holder holds any VP-Series ion exchange water purification cartridge. Connect to a pressurized water supply (up to 30 psi). Holds one cartridge (cartridges may be placed in sequence by mounting multiple holders). Quick, press-to-connect ends for ease of installation.

Made in the USA!



ORDERING GUIDE

SKU

UNIX-00

Includes:

- 6 feet of 1/4" OD Polyethylene tubing
- a shut-off valve
- inlet and outlet fittings
- mounting screws

Inlet and outlet connections must be installed by the user, and tubing cut to the desired length.

FILTERS & CONSUMABLES

SKU	Description
VP-17-4030	VP Series - High Purity Color-Changing DI Filter Cartridge (Std.)
VP-17-4040	VP Series - High Capacity Color-Changing DI Filter Cartridge (Std.)

Consumables sold seperately.